



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Tipack Rack DFL10ST, 0.1-10 µl, sterile, 10 x 96 pcs/pack

€136.00
plus VAT &
Shipping

Product Images



Description

Original sterile Gilson^{Diamond®} pipette tips with filter practically plugged in the refill rack. No more tedious plugging of tips. This saves time and is easy on the nerves.

DFL: Long tip version for volume 0.1–10 µl

TIPACK™, **Filter tips, sterile**

Gamma-sterilized, free of detectable RNases, DNases, DNA, RNA, ATP.

10 trays a 96 tips, individually packed with indicator label.

Additional Information

No.	CF-0423
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4403
EAN	4058072221942
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 µl
Volume MIN	0.1 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification with filter with box extra long RNase-free RNA-free DNase-free DNA-free ATP-free
Designation tip (Gilson)	DFL10ST
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Sterilization certificate from the manufacturer
Packaging type pipette tips (Gilson)	Sterilized TiPack™
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

