



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Tipack Rack DF30ST, 2-30 µl, sterile, 10 x 96 pcs/pack

**€123.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Original sterile Gilson<sup>Diamond®</sup> pipette tips with filter practically plugged in the refill rack. No more tedious plugging of tips. This saves time and is easy on the nerves.

### **TIPACK™, Filter tips, sterile**

Gamma-sterilized, free of detectable RNases, DNases, DNA, RNA, ATP.

10 trays a 96 tips, individually packed with indicator sticker.

## Additional Information

No.	CF-0428
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4404
EAN	4058072221997
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	30 µl
Volume MIN	2 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification with filter with box RNase-free RNA-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Designation tip (Gilson)	DF30ST
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Sterilization certificate from the manufacturer
Packaging type pipette tips (Gilson)	Sterilized TiPack™
for pipettes / dispensers of the manufacturers	Gilson

