



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex pipette tip Qualitix® 10 µL, with filter, thin graduated tip

**€170.77**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are not only suitable for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands

**Filter tips** include a hydrophobic 10 µm HDPE filter to block any contamination and effectively protect pipette and sample.

**Aerosol protection - no contamination !**

High tip purity as well as sterilization by gamma radiation are certified.

## Additional Information

No.	SR-0419
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7884
EAN	4058072369323
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	20 µl
Volume MIN	10 µl
sterile	Yes
Pipette tips properties	RNase-free RNA-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Model pipette tip	Socorex filter tips
for pipettes / dispensers of the manufacturers	Socorex

