



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix®, 100µL, with filter, slanted tip

€170.77
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are suitable not only for Socorex® pipettes but also for micropipettes of all brands. Different packaging units (bags, rack and the tipfill™ rack refill system), together with a wide range of tips with and without filters, offer enormous versatility.

Product advantages

- Extensive range, high material quality, autoclavable
- Bags, racks and tipfill™ rack refill system
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with micropipettes of all brands

Filter tips include a hydrophobic 10 µm HDPE filter to block any contamination and effectively protect pipette and sample.

Aerosol protection - no contamination !

High tip purity as well as sterilization by gamma radiation are certified.

Additional Information

No.	SR-0429
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8400
EAN	4058072369422
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	100 µl
Volume MIN	20 µl
sterile	Yes
Pipette tips properties	RNase-free RNA-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Model pipette tip	Socorex filter tips
for pipettes / dispensers of the manufacturers	Socorex

