



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Eppendorf epT.I.P.S.® Standard, Eppendorf Quality™, 0.2 – 5 mL L, 175 mm, purple, colorless tips, 300 tips (3 bags × 100 tips)**

**€40.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

With epT.I.P.S.® - the original Eppendorf "Totally Integrated Pipetting System" - the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. epT.I.P.S. Standard are high-quality original Eppendorf pipette tips, packed in resealable bags. They are available for volumes from 10 µL to 10 mL.

## Additional Information

No.	VB-0073
Article - USP	Original Eppendorf pipette tips, packed in resealable bags Available in sizes from 10 µL to 10 mL in purity grade Eppendorf Quality epT.I.P.S. Standard tips 200 µL, 300 µL and 1,000 µL are also available with the color codes yellow and
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
VGKL number	106930650
EAN	4058072176440
Gross weight	1.3
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	1.17000 kg
Color	purple
Material	Polypropylene (PP)
Volume MAX	5 ml
Volume MIN	0.2 ml
Country of origin	Germany
sterile	No
Packing width	0.34000 m
Packing height	0.18000 m
Packing depth	0.26000 m
Packaging volume	16000.000000 ccm
Customs tariff number	84799070
Length	175 mm
Degree of purity	Eppendorf Quality™
Packaging type Pipette tips (Eppendorf)	Standard

