



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**AHN Top-Line® Tip 1-200 µl, 50.8 mm, clear,
sterile, 10x96 pcs.**

€58.91 As low as
€45.31
**plus VAT &
Shipping**

Product Images



Description

The universal ones for all applications

Accurate pipetting without sample loss

Conical pipette outlet prevents drop formation at the tip opening

No risk of sample contamination

Pipette tips are made of autoclavable, non-wettable, pure polypropylene

Quick and easy volume control

Graduation lines on the pipette tip

- Free of RNase, DNase, pyrogens and PCR inhibitors
- Available in bags and rack boxes
- Securely seated and easy to discard pipette tips
- Optimal cone geometry with wide sealing zone
- Integrity of the product
- Sterile delivered product (beta irradiation)
- Compatible with neoPette, Eppendorf, Gilson, Brand, Biohit, Socorex, HTL

Additional Information

No.	AN-0033
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	AHN
old neoLab article no.	1-111-96-0
EAN	4058072195373
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	200 µl
Volume MIN	200 µl
sterile	Yes
Length	50.8 mm
Fluid behavior	hydrophobic
Packaging type pipette tips (AHN)	Rack

