




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



AHN myTip Gelloader 200 μ L, clear, rack box, carton/4x8x96

€129.95
plus VAT &
Shipping

Product Images

Description

Optimized fit and easy tip ejection - improved cone geometry for a secure fit Quick and easy volume control - graduation lines on the pipette tip Precise pipetting without sample loss - conical tip opening prevents drop formation at the tip outlet Use of high-quality material - tips made of autoclavable, non-wettable, pure polypropylene Easy identification of the pipette tip - colored racks allow quick and easy identification of the pipette tip as well as the appropriate pipette Product integrity - sterile version available (beta irradiation) Versions - available in bags and rack boxes Compatibility table on request

Additional Information

No.	AN-0041
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	AHN
EAN	4058072195397
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	3.70000 kg
Volume MAX	200 µl
Volume MIN	200 µl
Country of origin	Italy
sterile	No
Customs tariff number	84799070
Packaging type pipette tips (general)	Box

