



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand pipette tips BIO-CERT, 5-300 µl, neutral,
sterile, 10 x 96 pcs.**

**€77.10
plus VAT &
Shipping**

Product Images



Description

Brand Pipette tips for easy start.

Free of endotoxins, DNA, RNase and ATP.

Sterile pipette tips in certified BIO-CERT® quality.
With batch-related certificate.

Palletized (TipRack), 10 TipRacks à 96 tips

Refill unit, in environmentally friendly blister packaging. Supplied with transfer aid for transfer to an already autoclaved box.

Additional Information

No.	BN-1961
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Brand
old neoLab article no.	1-6527
EAN	4058072212001
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	300 µl
Volume MIN	5 µl
sterile	Yes
Pipette tips properties	with certification with box RNase-free Endotoxin-free DNA-free ATP-free
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Refill
Packing units tips	10 x 96
Model pipette tip	BRAND BIO-CERT Tips. Standard
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (BRAND)	TipRack
for pipettes / dispensers of the manufacturers	Biohit Brand Eppendorf Thermo Scientific neoLab

