



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand pipette tips BIO-CERT® palletized, 5-300  
µl, 10 x 96 pcs. + box, sterile**

**€80.90  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

BIO-CERT® pipette tips for easy entry.

Guaranteed sterile and free of endotoxins, DNA, RNase and ATP.

With batch-related certificate.

### TipBox sterile, 10 boxes à 96 tips

Box made of PP, with hinged hooded lid, stackable, multiple autoclavable at 121°C.

## Additional Information

---

No.	BN-1981
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Brand
old neoLab article no.	1-6522
EAN	4058072212186
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	300 µl
Volume MIN	5 µl
sterile	Yes
Pipette tips properties	with certification with box RNase-free Endotoxin-free DNA-free ATP-free
Fluid behavior	hydrophobic
Packing units tips	10 x 96
Model pipette tip	BRAND BIO-CERT Tips. Standard
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (BRAND)	TipBox
for pipettes / dispensers of the manufacturers	Biohit Brand Eppendorf Thermo Scientific

