



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand pipette tips BIO-CERT, 50-1000 µl, crystal,
sterile, 10 x 96 pcs.**

**€90.40
plus VAT &
Shipping**

Product Images



Description

Brand Pipette tips for easy start.

Free of endotoxins, DNA, RNase and ATP.

Sterile pipette tips in certified BIO-CERT® quality.

With batch-related certificate.

TipBox sterile, 10 boxes à 96 tips

Box made of PP, with hinged hooded lid, stackable, multiple autoclavable at 121°C.

Additional Information

| | |
|--|---|
| No. | BN-1982 |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Brand |
| old neoLab article no. | 1-6528 |
| EAN | 4058072212193 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Volume MAX | 1000 µl |
| Volume MIN | 50 µl |
| sterile | Yes |
| Pipette tips properties | with certification with box RNase-free Endotoxin-free DNA-free ATP-free |
| Fluid behavior | hydrophobic |
| Model pipette tip | BRAND BIO-CERT Tips. Standard |
| Type certificate | Batch certificate from the manufacturer |
| Packaging type pipette tips (BRAND) | TipBox |
| Packing units tips | 10 x 96 |
| for pipettes / dispensers of the manufacturers | Biohit Brand Eppendorf Thermo Scientific |

