



Thermo Scientific ClipTip™ pipette tips

Product Images

Description

ClipTip™ pipette tips Ensure a secure seal on every sample. Unique, innovative interlock technology ensures complete sealing of every channel with only minimal tip adhesion and ejection force. Description: ClipTip 1250, non-sterile Product line: ClipTip Volume (metric): 1250 µl Filtered: No Length (anglo): = 4.06 in. Compatibility: Thermo Scientific®; E1-ClipTip Compatibility: Thermo Scientific®; F1-ClipTip Format: Reload Tip Shape: Low Retention Packaging: 8 x 96 tips/insert (reload insert) Color Code: Turquoise Unit Size: Case of 768

Additional Information

No.	PO-0443
Manufacturer (Brand)	Thermo Scientific
EAN	4066292158160
Gross weight	1.34
Country of origin	Finnland
sterile	No
Packaging volume	0.02000000 ccm
Customs tariff number	84799070
Length	103 mm

