



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Pipette Tips, epTIPS Reloads, Quality™, 50-1250µl, 76mm, green

€86.20
plus VAT &
Shipping

Product Images



Description

epT.I.P.S.® Reloads, Eppendorf Quality™, 50 – 1,250 µL, 76 mm, green, colorless tips, 960 tips (10 trays × 96 tips) – Less waste compared to disposable racks – Reloads trays are autoclavable at up to 121 °C in the packaging and without foil for later use – Depending on tip size, trays are stacked against each other or stacked on top of each other in the refill system – Available in Eppendorf Quality and PCR clean purity grades – Color-coded trays for easy volume and pipette identification

Additional Information

| | |
|-----------------------|------------------|
| No. | VB-0169 |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Eppendorf |
| EAN | 4058072405014 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | green |
| Volume MAX | 1250 µl |
| Volume MIN | 50 µl |
| sterile | No |
| Length | 76 mm |
| Fluid behavior | hydrophobic |

