



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Eppendorf Pipette Tips, epT.I.P.S. Set 0.5–20µl

**€54.80**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

With epT.I.P.S.®-, the original Eppendorf "Totally Integrated Pipetting System"-, the tips are optimally matched to the Eppendorf pipettes and, of course, also meet the requirements of EN ISO 8655. This results in minimal attachment and ejection forces with maximum tightness. epT.I.P.S. pipette tips can also be used with pipettes from other manufacturers. Optimal wetting properties, high transparency and special, certified purity grades are the visible expression of our product and production philosophy. The epT.I.P.S. sets in Eppendorf Quality contain one box and 5 refill trays with tips. The refill system is completely autoclavable and available for volumes from 10 µL to 2.5 mL. Contamination-free transfer of trays into the stable work box System optimized for multichannel pipettes Color-coded trays for easy tip and pipette identification Tips can be placed directly from the refill trays The trays and boxes are autoclavable at 121 °C for later use epT.I.P.S. Box and epT.I.P.S. Set available in Eppendorf Quality

## Additional Information

No.	VB-0144
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072171322
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	20 µl
Volume MIN	0.5 µl
sterile	No
Fluid behavior	hydrophobic
Model pipette tip	Eppendorf epT.I.P.S.

