



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Pipette Tips, epTIPS Biopur, 0.1-5ml, 120mm, purple

€60.00
plus VAT &
Shipping

Product Images



Description

epT.I.P.S.® Racks, Biopur®, 0.1 – 5 mL, 120 mm, purple, colorless tips, 120 tips (5 racks × 24 tips). Replaces order no. 0030075200 – Eppendorf Biopur pipette tips in disposable racks have the highest biological purity level – Safe PCR clean, sterile, pyrogen-, ATP- and DNA-free – All pre-sterilized epT.I.P.S. and ep Dualfilter T.I.P.S. variants are available in racksVariants are available in racks – Sterile packaging complies with DIN EN ISO 11607 and DIN EN 868-2-10 standards – Each production batch of racks is continuously tested by an independent laboratory – Batch-specific certificates available at: www.eppendorf.com/certificates – Packaged in racks of 384, 96, 48 or 24 tips – For use in the pharmaceutical and food industries, molecular biology and cell technology – Modern, slim rack with resealable lid – Improved sustainability due to 35 % lower polypropylene consumption per rack

Additional Information

No.	VB-2614
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072425944
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	5 ml
Volume MIN	0.1 ml
sterile	Yes
Fluid behavior	hydrophobic
Packing units tips	5 x 24
Model pipette tip	Eppendorf epT.I.P.S.
Type certificate	Quality certificate
Packaging type Pipette tips (Eppendorf)	Racks
for pipettes / dispensers of the manufacturers	Eppendorf

