



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf epTIPS Standard 1-10 ml L, 243 mm
long, 2x100 pcs.**

**€85.80
plus VAT &
Shipping**

Product Images



Description

Top quality at an attractive price!

Original Eppendorf standard tips in resealable bags, autoclavable at 121°C.

epT.I.P.S.[®] Standard

Loosely packed in bag

Additional Information

No.	VB-0076
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6810
EAN	4058072173012
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 ml
Volume MIN	0.1 ml
sterile	No
Pipette tips properties	with certification graduated extra long
Fluid behavior	hydrophobic
Model pipette tip	Eppendorf epT.I.P.S.
Type certificate	Batch certificate from the manufacturer
Packaging type Pipette tips (Eppendorf)	Standard
Packaging type pipette tips (general)	Bag
Packing units tips	2 x 100
for pipettes / dispensers of the manufacturers	Eppendorf

