



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Eppendorf epTIPS Standard, 0.1-20 µl, 2x500 pcs.**

**€69.60**  
**plus VAT &**  
**Shipping**

## Product Images



## Description

---

Top quality at an attractive price!

Original Eppendorf standard tips in resealable bags, autoclavable at 121°C.

### **epT.I.P.S.<sup>®</sup> Standard**

Loosely packed in bag

## Additional Information

No.	VB-0078
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6806
EAN	4058072404352
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	20 µl
Volume MIN	0.1 µl
sterile	No
Pipette tips properties	with certification
Fluid behavior	hydrophobic
Model pipette tip	Eppendorf epT.I.P.S.
Type certificate	Batch certificate from the manufacturer
Packaging type Pipette tips (Eppendorf)	Standard
Packaging type pipette tips (general)	Bag
Packing units tips	2 x 500
for pipettes / dispensers of the manufacturers	Eppendorf

