



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epTIPS Standard, 2-200 µl, 2x500 pcs.

€42.90
plus VAT &
Shipping

Product Images



Description

Top quality at an attractive price!

Original Eppendorf standard tips in resealable bags, autoclavable at 121°C.

epT.I.P.S.[®] Standard

Loosely packed in bag

Additional Information

| | |
|------------------------------------------------|-----------------------------------------|
| No. | VB-0080 |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Eppendorf |
| old neoLab article no. | E-6802 |
| EAN | 4058072404376 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (yellow) |
| Volume MAX | 200 µl |
| Volume MIN | 2 µl |
| sterile | No |
| Pipette tips properties | with certification |
| Fluid behavior | hydrophobic |
| Model pipette tip | Eppendorf epT.I.P.S. |
| Type certificate | Batch certificate from the manufacturer |
| Packaging type Pipette tips (Eppendorf) | Standard |
| Packaging type pipette tips (general) | Bag |
| Packing units tips | 2 x 500 |
| for pipettes / dispensers of the manufacturers | Eppendorf |

