



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epTIPS Standard, 2-200 µl, 1000 pcs./bag

€42.90
plus VAT &
Shipping

Product Images



Description

Top quality at an attractive price!

Original Eppendorf standard tips in resealable bags, autoclavable at 121°C.

epT.I.P.S.[®] Standard

Loosely packed in a bag

Additional Information

No.	VB-0081
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6811
EAN	4058072174286
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	200 µl
Volume MIN	2 µl
sterile	No
Pipette tips properties	with certification
Fluid behavior	hydrophobic
Model pipette tip	Eppendorf epT.I.P.S.
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Eppendorf

