



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Eppendorf epTIPS Standard, 20–300 µl, 2x500
pcs.**

**€42.90
plus VAT &
Shipping**

Product Images



Description

Top quality at an attractive price!

Original Eppendorf standard tips in resealable bags, autoclavable at 121°C.

epT.I.P.S.[®] Standard

Loosely packed in bag

Additional Information

No.	VB-0082
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6807
EAN	4058072404383
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (yellow)
Volume MAX	300 µl
Volume MIN	20 µl
sterile	No
Pipette tips properties	with certification
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Bag
Packing units tips	2 x 500
Model pipette tip	Eppendorf epT.I.P.S.
Type certificate	Batch certificate from the manufacturer
Packaging type Pipette tips (Eppendorf)	Standard
for pipettes / dispensers of the manufacturers	Eppendorf

