



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Microman Tips CP10ST (sterile) , Tipack Rack 10 x 96 pcs/pack

€424.00
plus VAT &
Shipping

Product Images



Description

Microman® flasks with capillaries (CP) made of non-wetting PP, free of color additives.

Sterile Free of RNase, DNase, DNA, RNA, proteases.

TIPACK rack: 10 boxes of 96 tips.

Additional Information

No.	CF-0260
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4445
EAN	4058072220808
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Volume MAX	10 µl
Volume MIN	1 µl
sterile	Yes

