



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Microman tips CP25, Tipack Rack 10 x 96 pcs/pack

€405.00
plus VAT &
Shipping

Product Images



Description

Microman® flasks with capillaries (CP) made of non-wetting PP, free of color additives.

Non-sterile

TIPACKrack

10 boxes of 96 tips each

Additional Information

No.	CF-0252
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4378
EAN	4058072220723
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	25 µl
Volume MIN	3 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	CP25
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	1 x 960
for model pipette	Microman M25 E
Model pipette tip	Gilson MICROMAN Tips
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TiPack™
for pipettes / dispensers of the manufacturers	Gilson

