



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Microman tips CP250, Tipack Rack 10 x 96 pcs/pack

€405.00
plus VAT &
Shipping

Product Images



Description

Microman® flasks with capillaries (CP) made of non-wetting PP, free of color additives.

Non-sterile

TIPACKrack

10 boxes of 96 tips each

Additional Information

No.	CF-0254
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4369
EAN	4058072220747
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	250 µl
Volume MIN	50 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	CP250
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	1 x 960
for model pipette	Microman M250 E
Model pipette tip	Gilson MICROMAN Tips
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TiPack™
for pipettes / dispensers of the manufacturers	Gilson

