



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Microman tips CP250, 200 pcs/pack loose

€79.00
plus VAT &
Shipping

Product Images



Description

Microman® flasks with capillaries (CP) made of non-wetting PP, free of color additives.

Non-sterile

Loose **in bag**

200 tips

Additional Information

No.	CF-0257
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4372
EAN	4058072220778
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	250 µl
Volume MIN	50 µl
sterile	Yes
Pipette tips properties	with certification
Fluid behavior	hydrophobic
Designation tip (Gilson)	CP250
Packaging type pipette tips (general)	Tips and plungers in carton
Packing units tips	1 x 200
for model pipette	Microman M250 E
Model pipette tip	Gilson MICROMAN Tips
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Tips and plungers in carton
for pipettes / dispensers of the manufacturers	Gilson

