



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Microman Tips CP1000, Tipack Rack, 2 x 91 pcs/pack

€87.00
plus VAT &
Shipping

Product Images



Description

Microman® flasks with capillaries (CP) made of non-wetting PP, free of color additives.

Non-sterile

TIPACKrack

2 boxes of 91 tips each

Additional Information

| | |
|--|---|
| No. | CF-0285 |
| Typ Pipettenspitze | Direktverdrängerspitzen |
| Manufacturer (Brand) | Gilson |
| old neoLab article no. | 7-4373 |
| EAN | 4058072220877 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (milky) |
| Volume MAX | 1000 µl |
| Volume MIN | 100 µl |
| sterile | No |
| Pipette tips properties | with certification mit box |
| Designation tip (Gilson) | CP1000 |
| Fluid behavior | hydrophobic |
| Model pipette tip | Gilson MICROMAN Tips |
| Type certificate | Declaration of conformity (standardization) |
| Packaging type pipette tips (Gilson) | Autoclavable TiPack™ |
| Packaging type pipette tips (general) | Multiple racks/ trays each with box |
| Packing units tips | 1 x 182 |
| for model pipette | Microman M1000 E |
| for pipettes / dispensers of the manufacturers | Gilson |

