



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson TowerPack DL10, 0.1-10 µl, 10 x 96 pcs/pack

€56.00  
plus VAT &  
Shipping

### Product Images



## Description

---

### **Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.**

All pipette tips, Tipack Racks and Tower Racks are autoclavable at 121°C.

The color coding of the tip containers corresponds to the color coding of the individual Pipetman models, so that the right tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

### **TOWERPACK™**

Refill system for universal box, saves space and plastic, autoclavable

10 trays à 96 tips

## Additional Information

No.	CF-0396
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4352
EAN	4058072221775
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	20 µl
Volume MIN	0.1 µl
sterile	No
Pipette tips properties	with certification extra long
Designation tip (Gilson)	D10
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TowerPack™
Packaging type pipette tips (general)	Refill. stackable
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

