



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Gilson TowerPack D200, 2-200 µl, 10 x 96
pcs/pack**

**€47.00
plus VAT &
Shipping**



Description

Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.

All pipette tips, Tipack Racks and Tower Racks are autoclavable at 121°C.

The color coding of the tip containers corresponds to the color coding of the individual Pipetman models, so that the right tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

TOWERPACK™

Refill system for universal box, saves space and plastic, autoclavable

10 trays à 96 tips

Additional Information

No.	CF-0397
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4353
EAN	4058072221782
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	200 µl
Volume MIN	2 µl
sterile	Yes
Pipette tips properties	with certification
Designation tip (Gilson)	D200
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TowerPack™
Packaging type pipette tips (general)	Refill. stackable
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

