



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson TowerPack D1000, 200-1000 µl, 7 x 96 pcs/pack

€39.00
plus VAT &
Shipping

Product Images



Description

Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.

All pipette tips, Tipack Racks and Tower Racks are autoclavable at 121°C.

The color coding of the tip containers corresponds to the color coding of the individual Pipetman models, so that the right tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

TOWERPACK™

Refill system for universal box, saves space and plastic, autoclavable

Order No. 7-4355: 7 trays à 96 tips

Additional Information

No.	CF-0398
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neolab article no.	7-4355
EAN	4058072221799
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification RNase-free RNA-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Designation tip (Gilson)	D1000
Packaging type pipette tips (general)	Refill. stackable
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TowerPack™
for pipettes / dispensers of the manufacturers	Gilson

