



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Tipack Rack DL10, 0.1-20 µl, 10 x 96 pcs/pack

**€61.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**The refill variant for Gilson PIPETMAN® Tips Diamond® in blister packaging.**

Individually packaged racks with color coding compatible with Gilson pipette models.

For maximum cleanliness and easy handling.

Simply remove the cover foil of the packaging and transfer the ready-to-use rack with tips into the corresponding box.

Rack and tips are autoclavable; for this, the blister packaging must be removed before autoclaving

!

**PU:** 10 x 96 pcs.

## Additional Information

No.	CF-0443
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4492
EAN	4058072222147
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	20 µl
Volume MIN	0.1 µl
sterile	Yes
Pipette tips properties	with certification with box extra long
Designation tip (Gilson)	DL10
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Blister
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

