



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Blister Refill D200, 2-200 µl, 10 x 96 pcs/pack

€50.00
plus VAT &
Shipping

Product Images



Description

The refill variant for Gilson PIPETMAN® Tips Diamond® in blister packaging.

Individually packaged racks with color coding compatible with Gilson pipette models.

For maximum cleanliness and easy handling.

Simply remove the cover foil of the packaging and transfer the ready-to-use rack with tips into the corresponding box.

Rack and tips are autoclavable; for this, the blister packaging must be removed before autoclaving

!

PU: 10 x 96 pcs.

Additional Information

| | |
|------------------------------------------------|---------------------------------------------|
| No. | CF-0449 |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Gilson |
| old neoLab article no. | 7-4493 |
| EAN | 4058072222161 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (milky) |
| Volume MAX | 200 µl |
| Volume MIN | 2 µl |
| sterile | No |
| Pipette tips properties | with certification mit box |
| Fluid behavior | hydrophobic |
| Designation tip (Gilson) | D200 |
| Packaging type pipette tips (general) | Multiple racks/ trays each with box |
| Packing units tips | 10 x 96 |
| Model pipette tip | Gilson PIPETMAN Tips Diamond |
| Type certificate | Declaration of conformity (standardization) |
| Packaging type pipette tips (Gilson) | Blister |
| for pipettes / dispensers of the manufacturers | Gilson |

