



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Blister Refill D300, 20-300 µl, 10 x 96 pcs/pack

€50.00
plus VAT &
Shipping

Product Images



Description

The refill variant for Gilson PIPETMAN® Tips Diamond® in blister packaging.

Individually packed racks with color coding compatible with Gilson pipette models.

For maximum cleanliness and easy handling.

Simply remove the cover foil of the packaging and transfer the ready-to-use rack with tips into the appropriate box.

Rack and tips are autoclavable; therefore the blister packaging must be removed before autoclaving

!

PU: 10 x 96 pcs.

Additional Information

No.	CF-0455
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4494
EAN	4058072222192
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	300 µl
Volume MIN	20 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D300
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Blister
for pipettes / dispensers of the manufacturers	Gilson

