



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Gilson EasyPack D1200, 200–1200 µl, 5 x 200
pcs/pack**

**€31.00
plus VAT &
Shipping**

Product Images



Description

Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.

All pipette tips, Tipack Racks and Tower Racks are autoclavable at 121°C.

The color coding of the tip containers corresponds to the color coding of the individual Pipetman models, so that the right tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

EasyPack™

5 bags of 200 tips each

Additional Information

| | |
|--|---|
| No. | CF-0291 |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Gilson |
| old neoLab article no. | 7-4366 |
| EAN | 4058072220938 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (milky) |
| Volume MAX | 1200 µl |
| Volume MIN | 100 µl |
| sterile | No |
| Pipette tips properties | with certification |
| Fluid behavior | hydrophobic |
| Designation tip (Gilson) | D1200 |
| Packaging type pipette tips (general) | Bag |
| Packing units tips | 5 x 200 |
| Model pipette tip | Gilson PIPETMAN Tips Diamond |
| Type certificate | Declaration of conformity (standardization) |
| Packaging type pipette tips (Gilson) | EasyPack™ |
| for pipettes / dispensers of the manufacturers | Gilson |

