



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson EasyPack D10ml, 1-10 ml, 4 x 50 pcs/pack

**€55.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.**

All pipette tips, Tipack Racks and Tower Racks are autoclavable at 121°C.

The color coding of the tip containers corresponds to the color coding of the individual Pipetman models, so that the right tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

### **EasyPack™**

Order No. 7-4341: 4 bags of 50 tips each

## Additional Information

---

No.	CF-0295
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4341
EAN	4058072220945
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 µl
Volume MIN	1 µl
sterile	No
Pipette tips properties	with certification
Fluid behavior	hydrophobic
Designation tip (Gilson)	D10ml
Packaging type pipette tips (general)	Bag
Packing units tips	4 x 50
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	EasyPack™
for pipettes / dispensers of the manufacturers	Gilson

