



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson EcoPack D10, 0.1-10 µl, 10 x 1000 pcs/pack

€281.00
plus VAT &
Shipping

Product Images



Description

Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.

All pipette tips, Tipack Racks and Tower Racks are autoclavable at 121°C.

The color coding of the tip containers corresponds to the color coding of the individual Pipetman models, so that the right tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

EcoPack™

Cost-effective bulk pack

10 bags of 1000 tips each

Additional Information

No.	CF-0330
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4334
EAN	4058072221171
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 µl
Volume MIN	0.1 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D10
Packaging type pipette tips (general)	Bag
Packing units tips	10 x 1000
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	EcoPack™
for pipettes / dispensers of the manufacturers	Gilson

