



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Reload Pack D200, 2-200 µl, 10 x 96 pcs/pack

€39.00
plus VAT &
Shipping

Product Images



Description

Original Gilson Diamond® pipette tips for all Gilson microliter pipettes in the Reload Pack. The Reload Pack is the easiest and fastest way to fill an empty TiPack box. Environmentally friendly because 70% less plastic packaging. The tips are made of transparent polypropylene. Racks are color coded for better differentiation.

Rack and tips are autoclavable; for this, the reload packaging must be removed before autoclaving !

PU: 10 x 96 pcs.

Additional Information

No.	CF-0392
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4408
EAN	4058072188368
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	200 µl
Volume MIN	2 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D200
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Reload Pack
for pipettes / dispensers of the manufacturers	Gilson

