



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Tipack Rack D10ml, 1-10 ml, 8 x 30 pcs/pack

**€140.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Original Gilson<sup>Diamond®</sup> pipette tips practically plugged in the rack. No more tedious plugging of tips. This saves time and is easy on the nerves.

All tips are made of sturdy polypropylene and are autoclavable. Available in different volumes.

TIPACK™

Autoclavable, color-coded sealed.

8 trays of 30 tips each

## Additional Information

No.	CF-0297
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4401
EAN	4058072220969
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 µl
Volume MIN	1 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D10ml
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	8 x 30
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Sterilized TiPack™
for pipettes / dispensers of the manufacturers	Gilson

