



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Tipack Rack D200, 2-200 µl, 10 x 96 pcs/pack

€54.00
plus VAT &
Shipping

Product Images



Description

Original Gilson^{Diamond®} pipette tips practically plugged in the rack. No more tedious plugging of tips. This saves time and is easy on the nerves.

All tips are made of sturdy polypropylene and are autoclavable. Available in different volumes.

TIPACK™

Autoclavable, color-coded sealed.

10 trays à 96 tips

Additional Information

No.	CF-0425
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4316
EAN	4058072221966
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	200 µl
Volume MIN	2 µl
sterile	Yes
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D200
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TiPack™
for pipettes / dispensers of the manufacturers	Gilson

