



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Tipack Rack D300, 20-300 µl, 10 x 96 pcs/pack

€54.00
plus VAT &
Shipping

Product Images



Description

Original Gilson^{Diamond®} pipette tips practically plugged in the rack. No more tedious plugging of tips. This saves time and is easy on the nerves.

All tips are made of sturdy polypropylene and are autoclavable. Available in different volumes.

The tips with a volume of 0.1-10 µl are also available in extra long.

TIPACK™

Autoclavable, color-coded sealed.

10 trays à 96 tips

Additional Information

| | |
|--|---|
| No. | CF-0429 |
| Typ Pipettenspitze | Pipettenspitzen |
| Manufacturer (Brand) | Gilson |
| old neoLab article no. | 7-4338 |
| EAN | 4058072222000 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (milky) |
| Volume MAX | 300 µl |
| Volume MIN | 20 µl |
| sterile | Yes |
| Pipette tips properties | with certification mit box |
| Fluid behavior | hydrophobic |
| Designation tip (Gilson) | D300 |
| Packaging type pipette tips (general) | Multiple racks/ trays each with box |
| Packing units tips | 10 x 96 |
| Model pipette tip | Gilson PIPETMAN Tips Diamond |
| Type certificate | Declaration of conformity (standardization) |
| Packaging type pipette tips (Gilson) | Autoclavable TiPack™ |
| for pipettes / dispensers of the manufacturers | Gilson |

