



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Tipack Rack D1200, 100-1200 μ l, 10 x 96 pcs/pack

€57.00
plus VAT &
Shipping

Product Images



Description

Original Gilson^{Diamond®} pipette tips practically plugged in the rack. No more tedious plugging of tips. This saves time and is easy on the nerves.

All tips are made of sturdy polypropylene and are autoclavable. Available in different volumes.

TIPACK™

Autoclavable, color-coded sealed.

10 trays à 96 tips

Additional Information

No.	CF-0438
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4312
EAN	4058072222093
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1200 µl
Volume MIN	100 µl
sterile	No
Pipette tips properties	with certification mit box
Designation tip (Gilson)	D1200
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TiPack™
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

