



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson TowerPack D10ST, 0.1-10 µl, 10 x 96 pcs/pack

**€62.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

When it has to be more again: **Original Gilson Diamond® pipette tips, sterile**, for all Gilson microliter pipettes in **Tower Pack**. No more tedious plugging. This saves time and is easy on the nerves. In addition, the tower rack saves up to 60% space.

Tips made of transparent polypropylene, each tip marked in two places according to its volume.

### **TOWERPACK™, sterile**

Refill system for universal box

10 trays of 96 tips each

Gamma sterilized, free of detectable RNases, DNases, DNA, RNA, ATP.

All trays in one package, not individually wrapped.

## Additional Information

|                                                |                                                                                  |
|------------------------------------------------|----------------------------------------------------------------------------------|
| No.                                            | CF-0400                                                                          |
| Typ Pipettenspitze                             | Pipettenspitzen                                                                  |
| Manufacturer (Brand)                           | Gilson                                                                           |
| old neolab article no.                         | 7-4348                                                                           |
| EAN                                            | 4058072221812                                                                    |
| Transport temperature                          | Room temperature                                                                 |
| autoclavable at 121°C                          | Yes                                                                              |
| Color                                          | transparent (milky)                                                              |
| Volume MAX                                     | 10 µl                                                                            |
| Volume MIN                                     | 0.1 µl                                                                           |
| sterile                                        | No                                                                               |
| Pipette tips properties                        | pyrogen-free with certification RNase-free RNA-free DNase-free DNA-free ATP-free |
| Designation tip (Gilson)                       | D10ST                                                                            |
| Fluid behavior                                 | hydrophobic                                                                      |
| Model pipette tip                              | Gilson PIPETMAN Tips Diamond                                                     |
| Type certificate                               | Declaration of conformity (standardization)                                      |
| Packaging type pipette tips (Gilson)           | Autoclavable TowerPack™                                                          |
| Packaging type pipette tips (general)          | Refill. stackable                                                                |
| Packing units tips                             | 10 x 96                                                                          |
| for pipettes / dispensers of the manufacturers | Gilson                                                                           |

