



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson TowerPack DL10ST, 0.1-20 µl, 10 x 96 pcs/pack

€62.00
plus VAT &
Shipping

Product Images



Description

When it has to be more again: **Original Gilson Diamond® pipette tips, sterile**, for all Gilson microliter pipettes in **Tower Pack**. No more tedious plugging. This saves time and is easy on the nerves. In addition, the tower rack saves up to 60% space.

Tips made of transparent polypropylene, each tip marked in two places according to its volume.

DL: Long tip version for volume 0.1–20 µl

TOWERPACK™, sterile

Refill system for universal box

10 trays of 96 tips each

Gamma sterilized, free of detectable RNases, DNases, DNA, RNA, ATP.

All trays in one package, not individually wrapped.

Additional Information

No.	CF-0401
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4349
EAN	4058072171780
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	20 µl
Volume MIN	0.1 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification extra long RNase-free RNA-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Designation tip (Gilson)	DL10ST
Packaging type pipette tips (general)	Refill. stackable
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TowerPack™
for pipettes / dispensers of the manufacturers	Gilson

