



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Gilson TowerPack D1000ST, 100-1000  $\mu$ l, 7 x 96  
pcs/pack**

**€47.00  
plus VAT &  
Shipping**



## Description

---

When it has to be more again: **Original Gilson Diamond® pipette tips, sterile**, for all Gilson microliter pipettes in **Tower Pack**. No more tedious plugging. This saves time and is easy on the nerves. In addition, the tower rack saves up to 60% space.

Tips made of transparent polypropylene, each tip marked in two places according to its volume.

### **TOWERPACK™, sterile**

Refill system for universal box

7 trays of 96 tips each

Gamma sterilized, free of detectable RNases, DNases, DNA, RNA, ATP.

All trays in one package, not individually wrapped.

## Additional Information

No.	CF-0403
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neolab article no.	7-4358
EAN	4058072221836
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification RNase-free RNA-free DNase-free DNA-free ATP-free
Designation tip (Gilson)	D1000ST
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Autoclavable TowerPack™
Packaging type pipette tips (general)	Refill. stackable
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

