



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Blister Refill DL10ST, 0.1-20 µl, 10 x 96 pcs/pack

€68.00
plus VAT &
Shipping

Product Images



Description

Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.

The color coding of the racks corresponds to the color coding of the individual Pipetman models, so that the matching tip can be found in no time.

Tips made of crystal clear polypropylene, each tip marked in two places according to its volume.

Blister Refill™

Sterile, color-coded sealed.

10 trays individually packed with 96 tips each

Additional Information

No.	CF-0444
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4496
EAN	4058072222154
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	20 µl
Volume MIN	0.1 µl
sterile	No
Pipette tips properties	pyrogen-free with certification with box extra long RNase-free RNA-free DNase-free DNA-free
Designation tip (Gilson)	DFLI0ST
Fluid behavior	hydrophobic
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Blister
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
for pipettes / dispensers of the manufacturers	Gilson

