




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Blister Refill D300ST, 20-300 µl, 10 x 96 pcs/pack

**€58.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **Original Gilson pipette tips for all Gilson PIPETMAN microliter pipettes.**

The color coding of the racks corresponds to the color coding of the individual pipetman models, so that the matching tip can be found in no time.

Tips made of crystal-clear polypropylene, each tip marked in two places according to its volume.

#### **Blister Refill™**

Sterile, color-coded sealed.

10 trays individually packed with 96 tips each

## Additional Information

No.	CF-0456
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4498
EAN	4058072222208
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	300 µl
Volume MIN	20 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification with box RNase-free RNA-free DNase-free DNA-free
Fluid behavior	hydrophobic
Designation tip (Gilson)	DF300ST
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Blister
for pipettes / dispensers of the manufacturers	Gilson

