



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson Tipack Rack D10ST, 0.1-10 µl, 10 x 96 pcs/pack

€73.00
plus VAT &
Shipping

Product Images



Description

Original Gilson^{Diamond®} pipette tips, sterile, conveniently rack-mounted.

Gamma-sterilized; detection by indicator sticker.

Free of detectable RNases, DNases, DNA, RNA, ATP.

No more tedious plugging of tips. Saves time and is easy on the nerves.

TIPACK™

Sterile, color-coded sealed.

10 trays à 96 tips

Additional Information

No.	CF-0417
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4313
EAN	4058072172152
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 µl
Volume MIN	0.1 µl
sterile	No
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	DF10ST
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Sterilized TiPack™
for pipettes / dispensers of the manufacturers	Gilson

