



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Tipack Rack D1000ST, 100-1000 µl, 10 x 96 pcs/pack

**€63.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Original Gilson<sup>Diamond®</sup> pipette tips, sterile, conveniently rack-mounted.

Gamma-sterilized; detection by indicator sticker.

Free of detectable RNases, DNases, DNA, RNA, ATP.

No more tedious plugging of tips. Saves time and is easy on the nerves.

TIPACK™

Sterile, color-coded sealed.

10 trays à 96 tips

## Additional Information

No.	CF-0434
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4345
EAN	4058072222055
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	Yes
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D1000ST
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Sterilized TiPack™
for pipettes / dispensers of the manufacturers	Gilson

