



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Tipack Rack D1200ST, 100-1200 µl, 10 x 96 pcs/pack

**€63.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Original Gilson<sup>Diamond®</sup> pipette tips, sterile, conveniently rack-mounted.

Gamma-sterilized; detection by indicator sticker.

Free of detectable RNases, DNases, DNA, RNA, ATP.

No more tedious plugging of tips. Saves time and is easy on the nerves.

TIPACK™

Sterile, color-coded sealed.

10 trays à 96 tips

## Additional Information

No.	CF-0439
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4346
EAN	4058072222109
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	1200 µl
Volume MIN	100 µl
sterile	Yes
Pipette tips properties	with certification mit box
Fluid behavior	hydrophobic
Designation tip (Gilson)	D1200ST
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Sterilized TiPack™
for pipettes / dispensers of the manufacturers	Gilson

