



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® MicroFlex microcapillary pipette tips 0.5–200 µl, sterile, 200 pcs/pack

€64.40
plus VAT &
Shipping

Product Images



Description

Microcapillary pipette tips made of PP, universally applicable with microliter pipettes from Eppendorf, Brand, Gilson, Finnpiquette and Socorex.

Especially for gel electrophoresis.

- Round tips with 0.57 mm outer Ø, for 0.75, 1.0 and 1.5 mm gels.

RNase /DNase-free

Color: transparent

Tips not autoclavable.

Additional Information

No.	1-6401
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	neoLab
EAN	4058072071639
autoclavable at 121°C	Yes
diameter	0,57 mm
Color	transparent (crystal clear)
Volume MAX	200 µl
Volume MIN	0.5 µl
sterile	Yes
Pipette tips properties	RNase-free DNA-free
Fluid behavior	hydrophobic

