



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® pipette tips yellow, 1-200 µl, 1000
pcs./bag**

**€13.50
plus VAT &
Shipping**

Product Images



Description

High-quality disposable tips made of PP in common sizes, in a bag.
Suitable for most microliter pipettes. (Eppendorf, Brand, Gilson)

Additional Information

No.	E-1000
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	neoLab
EAN	4058072062095
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (orange)
Volume MAX	200 µl
Volume MIN	1 µl
sterile	No
Packaging type pipette tips (general)	Bag
Packing units tips	1 x 1000
for pipettes / dispensers of the manufacturers	Eppendorf Brand Gilson

