



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand pipette tips, 2 - 200 µl, PP, colorless, in bag, length 50 mm, non-sterile

Special Price

€11.65 was

€23.30

**plus VAT &
Shipping**

Product Images



Description

Pipette tips and filter tips from BRAND offer reliable quality for precise and demanding analyses. The raw materials are free from the additives DiHEMA and oleamide. Tips are manufactured and automatically packaged at BRAND under controlled clean room conditions. No release agents or demolding agents are used. Palletized tips and filter tips up to 1000 µl are free of DNA (<40 fg), RNase (<8.6 fg), endotoxins (<1 pg) and ATP (<1 fg). The BRAND tip range includes pipette tips for volumes from 0.1 µl to 10 ml and filter tips for volumes from 0.1 µl to 1000 µl. They are tested for Transferpette® microliter pipettes from BRAND and for many other pipettes from other manufacturers. For more information, please refer to the Selection Guide. BRAND pipette tips and filter tips are autoclavable at 121 °C (2 bar) according to DIN EN 285 and CE-marked according to IVD Directive 98/79 EC. Sterile pipette and filter tips are manufactured in certified BIO-CERT® quality. They are sterilized with β-rays according to ISO 11 137 and AAMI guidelines.

Additional Information

No.	BN-1942
Typ Pipettenspitze	Pipettenspitzen
Article - USP	Cleanroom quality for contamination-free work. Perfect fit on the microliter pipette Transferpette® and many other pipettes. The top system for all requirements.
Manufacturer (Brand)	Brand
VGKL number	104932008
EAN	4058072178666
Gross weight	0.372
HSNumber	39269097
Transport temperature	Room temperature
Storage temperature	Room temperature
autoclavable at 121°C	Yes
Net weight	0.27800 kg
Volume MAX	200 µl
Volume MIN	2 µl
sterile	No
Country of origin	Germany
Packing width	0.13000 m
Packing height	0.11000 m
Packing depth	0.29000 m
Packaging volume	4147.000000 ccm
Customs tariff number	39269097
Packaging type pipette tips (BRAND)	Bag

