



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex pipette tip Qualitix® for gel, 200µL

€113.88
plus VAT &
Shipping

Product Images



Description

Socorex Qualitix® pipette tips meet the highest quality standards, offer high flexibility and reliable compatibility. The pipette tips are not only suitable for Socorex® pipettes but also for micropipettes of all brands.

Product advantages

- Comprehensive range, high material quality, autoclavable
- Guaranteed consistency, batch after batch
- Certified nuclease-, endotoxin- and ATP-free
- Reduced packaging requires less storage space in the lab
- Versatile compatibility with all brands of micropipettes

Capillary tips for loading gels feature a long, soft end.

Additional Information

No.	SR-0440
Typ Pipettenspitze	Pipettenspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-8425
EAN	4058072369538
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	200 µl
Volume MIN	200 µl
sterile	No
Length	82.5 mm
Pipette tips properties	autoclavable
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Type certificate	Batch certificate from the manufacturer

